

BEAUTY AND YOUTH

50 PROGRAMS

SYNERGY OF QUANTUM AND
ELECTROMAGNETIC
TECHNOLOGY

POSSIBILITY OF USE
WITHOUT A QUANTUM NOZZLE



**Deta Elis
Egypt**
Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

BEAUTY AND YOUTH



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com

DEVITA COSMO FUNCTIONS

Takes care of the beauty of skin, hair and nails

Regulates the functioning of organs and systems

Improves the appearance of the skin

Restores the endocrine system
of women and men

Starts the rejuvenation process

Eliminates inflammation and pigmentation spots

Helps smooth out wrinkles

Restores skin tone and elasticity



**Deta Elis
Egypt**

Bio-Holistic Balance

2 Khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com

DEVITA COSMO FUNCTIONS



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com

DEVITA COSMO NEW - LIST OF PROGRAMS

1. Your protection
2. Regulation of blood circulation
3. Regulation of digestion
4. Cleansing the body
5. Live peacefully
6. Regulation of endocrine
7. Cleanse your kidneys
8. Cleanse the lymph
9. Clear gaze Female youth
10. Regulation of the heart
11. Stop colds
12. Easy breathing
13. Sleep regulation
14. Spiritual harmony
15. Stop joint inflammation
16. Healthy metabolism
17. Skin care
18. Beautiful nails and hair
19. The kidneys are fine
20. Vegetative is normal
21. Calm your nerves
22. Strong memory
23. Ears without pain
24. Healthy throat
25. Normalize the thyroid gland
26. Regulation of the stomach
27. Healthy muscles
28. Regeneration of injuries
29. Normalize the intestines
30. My head doesn't hurt!
31. Help for the bronchi
32. Male youth
33. Female youth
34. Like clockwork
35. Key of Life
36. Control your sugar levels
37. Stop pain (joints)
38. Strong joints
39. Tissue respiration
40. Prostate health
41. Concentration
42. Speed reading
43. I'm learning to count
44. Spleen support
45. Smoothing wrinkles
46. Facial tone
47. Vision correction
48. Hair beauty
49. Lightness of legs
50. Manicurist



**Deta Elis
Egypt**
Bio-Holistic Balance

2 Khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

DEVITA COSMO NEW

- 1.
- 2.
- 3.
- 4.
- 5.
- 6.
- 7.
- 8.
- 9.
- 10.
- 11.
- 12.
- 13.
- 14.
- 15.
- 16.
- 17.
- 18.
- 19.
- 20.
- 21.

- 22.
- 23.
- 24.
- 25.
- 26.
- 27.
- 28.
- 29.
- 30.
- 31.
- 32.
- 33.
- 34.
- 35.
- 36.
- 37.
- 38.
- 39.
- 40.
- 41.
- 42.

- 43.
- 44.
- 45.
- 46.
- 47.
- 48.
- 49.
- 50.



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com



THE DEVICE IS ATTACHED WITH DOCUMENTATION DESCRIPTION OF THE PROGRAMS, METHOD AND SCOPE OF THEIR APPLICATION USING A QUANTUM NOZZLE
6) OPERATING PRINCIPLE OF THE QUANTUM NOZZLE



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com



OPERATING PRINCIPLE OF THE QUANTUM NOZZLE

QUANTUM THERAPY -- is the effect of light and color on biological

Important points in the human body. Light, as a physical factor, consists of a large number of waves of different frequencies, so each color has its own frequency.

Each frequency has its own effect on the human body, which is why, to fully treat any condition, it's necessary to apply several frequencies and colors of light to specific areas of the skin. These areas are targeted with light waves that act as regulators of energy flow within the body.



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

OPERATING PRINCIPLE OF THE QUANTUM NOZZLE



Deta Elis
Egypt
Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

OPERATING PRINCIPLE OF THE QUANTUM NOZZLE

The key parameter is the radiation wavelength. Therefore, inexpensive LEDs have been developed for therapeutic purposes, producing red (660 + 15 nm wavelength), green (540 nm wavelength), blue (470 nm wavelength), and infrared (840-950 nm wavelength) radiation. The observed effect is similar to laser therapy.

Developers of equipment and application methods emphasize the technological advantages of this direction:

- significant reduction in the doses of prescribed pharmaceuticals;
- reduction of the rehabilitation period by 2-3 times;
- no side effects, non-invasive and environmentally safe.



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

OPERATING PRINCIPLE OF THE QUANTUM NOZZLE



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com




**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

PRIMARY SPECTRAL COLORS





**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,


Masr El Gedida Heliopolis

T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com



ПРИНЦИП ДЕЙСТВИЯ КВАНТОВОЙ НАСАДКИ

Рекомендуемые области применения аппарата в зависимости от длины волны и мощности излучения: длина волны 670 нм, мощность излучения 60 мВт (красный свет) применяется для расширения сосудов и восстановления кровотока; длина волны 505 нм, мощность излучения 12 мВт (зеленый свет) применяется при дерматологических проблемах; длина волны 440 нм, мощность излучения 12 мВт (синий свет), применяется при наличии гнойных и инфицированных ран.

Эффект от работы электромагнитного устройства проявляется после 3-5 процедур с экспозицией пять минут на частоту. При хронических проблемах эффект наблюдается после 10-20 процедур с ежедневным сеансом длительностью 10-30 мин.



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302

info4u@detaelisegypt.com
www.detaelisegypt.com

OPERATING PRINCIPLE OF THE QUANTUM NOZZLE

Recommended areas of application of the device

Depending on the wavelength and radiation power: wavelength 670 nm, radiation power 60 mW (red light) is used to dilate blood vessels and restore blood flow; wavelength 505 nm, radiation power 12 mW (green light) is used for dermatological problems; wavelength 440 nm, radiation power 12 mW (blue light) is used in the presence of purulent and infected wounds.

The effect of the electromagnetic device is evident after 3-5 treatments with a five-minute exposure per frequency. For chronic problems, the effect is observed after 10-20 treatments with daily sessions lasting 10-30 minutes.



**Deta Elis
Egypt**
Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

OPERATING PRINCIPLE OF THE QUANTUM NOZZLE



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
info4u@detaelisegypt.com
www.detaelisegypt.com

BEAUTY AND YOUTH



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302
infoAu@detaelisegypt.com
www.detaelisegypt.com

A unique device for self-administered cosmetic procedures. A fundamentally new approach to solving skin and appearance problems. It ensures high efficiency through the simultaneous reflex action of light quanta on nerve endings and the electromagnetic healing effect on the body's organs and systems.

The effect of beauty and rejuvenation is based on the principle of bioresonance and the reverse-reflex effect of quantum radiation.

EFFECTIVE. SAFE. PHYSIOLOGICAL.

EASY TO USE.



- **New Ritm+ power level**
- **New wireless Bluetooth nozzle**
- **Modern power system**
- **A new standard in comfort and efficiency**



**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt

Tell: +227747259 & +201023633302

info4u@detaelisegypt.com

www.detaelisegypt.com





**Deta Elis
Egypt**

Bio-Holistic Balance

2 khalil Al Arousi,
Masr El Gedida Heliopolis
T.K 11757 Cairo – Egypt
Tell: +227747259 & +201023633302

info4u@detaelisegypt.com
www.detaelisegypt.com



